


# GL SPECTROSOFT AUTOMATION

## add-on for GL SPECTROSOFT

It is an intelligent tool which communicates with and manages peripheral devices connected to your measurement system.



### GL SPECTROSOFT AUTOMATION

Instead of the need for setting up parameters in each plugged in device separately, this smart additional software cooperates with and runs all of them. Thanks to this, in one program on your computer, you can set sequences of operations such as: cooling temperature, current and voltage before a measurement. GL SPECTROSOFT AUTOMATION can collect data from external devices and link them together with measurement results.

GL SPECTROSOFT AUTOMATION automates complex operations, supplements and accelerates measurements of additional information (eg. the luminous efficacy, temperature, etc.) and minimises errors caused by misalignment of the equipment. It is a novelty on the market and a revolutionary solution in light measurement.

This is perfect tool for laboratories, where there is the need for complex operations. The AUTOMATION ADD-ON supports programmable: power supplies, current sources, multimeters and temperature controllers. Automation software can communicate with peripheral devices using a custom-made plug-in.

GL SPECTROSOFT AUTOMATION is available with GL SPECTROSOFT PRO or LAB versions only. Integration of new devices is available upon request.

