### **GL GONIOPHOTOMETER**

GLG A 50-1800

The GLG A 50-1800 is Type A goniometer system compliant with the CIE 121-1996 and IESNA LM-75-01 standards regulating far-field photometric and colorimetric measurement systems. It is designed for the testing of light sources specific to automotive, railway or similar industries basing on regulations specified by UN/ECE and SAE/FMVSS.



- H, V axis and X, Y, Z direction movement with servo motors
- Angular resolution of 0.002°
- Universal mounting table with position adjustability in three axes
- 10 inch LCD touch screen and digital manipulator
- Compatible with photometer, spectroradiometer, colorimeter and retroreflectometer
- Reduced weight for easier transport and installation
- Designed for safety and convenience of operation



# **GL GONIOPHOTOMETER**

### GLG A 50-1800

#### APPLICATION

Compliance testing of lamps and signals to UN/ECE and SAE /FMVSS Regulations.

This is a Type A system compliant with the CIE 121-1996 and IESNA LM-75-01 standards for far-field photometric and colorimetric measurement systems.

TECHNICAL DATA SHEET	
CIE Goniometer type	<ul> <li>Far Field Type A with H, V axis and x,y,z direction movement</li> <li>DUT positioning in X,Y axes on a motorized table</li> <li>5 axis servo motors with absolute encoders</li> </ul>
H axis movement	Angular range ± 180°. Nominal torque 98 Nm. Speed up to 50°/s
H axis resolution	0.002°
H axis reproducibility	0.05° * (*at nominal load)
V axis movement	Angular range ± 100°. Nominal torque 555 Nm. Speed up to 10°/s
V axis resolution	0.002°
V axis reproducibility	0.05° * (*at nominal load)
Z axis movement	Linear range 0-600 mm. Lifting capability up to 1500 N
Z axis reproducibility	70 μm
X axis movement	Linear range ± 300 mm. Speed up to 40 mm/s
Y axis movement	Linear range ± 150 mm. Speed up to 40 mm/s
X,Y axis reproducibility	70 μm
DUT mounting table	Square 500x500 mm. T-slot array 5x5 in the Bosch Rexroth system
Maximum DUT dimensions	≤ 1800 mm (symmetrically positioned)
Goniometer dimensions	2120(L) x 1780(H) x 800(W) mm
Optical axis height	1500 mm
Minimum room height	2600 mm
Minimum room width	4000 mm
Maximum load	50 kg
Goniometer weight	600 kg
Power supply	AC 3-phase 110-230V, 3000 W
Controllers	<ul> <li>Internal controller for all motors</li> <li>PC connectivity through LAN network</li> <li>Remote manufacturer support capability</li> <li>10 inch LCD touch display on the external control panel</li> <li>Manual all axis controller for operations during DUT fixing</li> </ul>
Safety	External control panel with emergency stop button     Laser scanner
Sensor type	<ul> <li>GL Photometer 3.0 LS + Flicker</li> <li>GL SPECTIS 1.3 LS</li> <li>GL SPECTIS 4.0</li> </ul>
Product no.	202522

**Note:** Instrument, firmware and software specification are subject to change without prior notice. All information included in GL OPTIC datasheets and product information available in any form are carefully prepared and included information believed to be true. Please note that discrepancies may occur due to text and/or other errors or changes in the available technology. We advise to contact GL Optic before the use of the product to obtain the latest product specification. Published on March 13th 2023.







# TECHNICAL SHEET





